



Recombinant Arabidopsis thaliana F-box protein SKIP3 (SKIP3)

Product Code	CSB-EP872324DOA
Abbreviation	SKIP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FV02
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MIRLRTLIEK NRSESRKKIV SSTLDSLPEG CISNIISFTS PEDACVAAAV SKIFESAVKS DIVWEKFLPT DYESLITPSR VFSSKKELYF SLCNDPLLIE DGKMSLWLEK ASGKRCIMLS ATAMNLSSMA DMSQRFLWIP CPESRFETVA ALREAYRFEF NCRMNTRVLS LRTRYSVYIV FKKADNWCGF KGVSIEAVVG IVGEESFRSF ICFDTHGKGQ ARKRKVAKP ELREDGWMET EIGEFYNEGG LMSSDEVEIS TVEGKYAQQK RGLVILGIEI RPAKILEEFS VCRL
Source	E.coli
Target Names	SKIP3
Protein Names	Recommended name: F-box protein SKIP3 Alternative name(s): F-box protein PP2-B9 Protein PHLOEM PROTEIN 2-LIKE B9 SKP1-interacting partner 3
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.