



Recombinant Arabidopsis thaliana Transcription factor TCP1 (TCP1)

Product Code	CSB-EP872344DOA
Abbreviation	TCP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FYG7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSSSTNDYND GNNNGVYPLS LYLSSLSGHQ DIIHNPYNHQ LKASPGHMVS AVPESLIDYM AFKSNNVVNQ QGFEPPEVSK EIKKVVKKDR HSKIQTAQGI RDRRVRLSIG IARQFFDLQD MLGFDKASKT LDWLLKKSrk AIKEVVQAKN LNNDDEDFGN IGGDVEQEEE KEEDDNGDKS FVYGLSPGYG EEEVVCEATK AGIRKKKSEL RNISSKGLGA KARGKAKERT KEMMAYDNPE TASDITQSEI MDPFKRSIVF NEGEDMTHLF YKEPIEEFDN QESILTNTMTL PTKMGQSYNQ NNGILMLVDQ SSSSNYNTFL PQNLDSYDQ NPFHDQTLVY VTDKNFPKKGK VWIQDSFVN
Source	E.coli
Target Names	TCP1
Protein Names	Recommended name: Transcription factor TCP1
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.