



Recombinant Arabidopsis thaliana F-box/FBD/LRR-repeat protein At1g13570 (At1g13570)

Product Code	CSB-MP872349DOA
Abbreviation	At1g13570
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FZ70
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MGESPFDISD LPQSIENIL TRLSIRDAIR TSVLSSKWRY KWSTLTDLVF DEKCVSPSND RCVVETNLVR FITGVLLHQQ GPIHKFQLST SFKQCRPDID QWLLFLSRNG IKELVLKLG EFRVPAFLF NCLKLTCEL CHCEFDPPQY FKGFSYLKSL NLHQILVAPE VIESLISGCP LLEFLSLSYF DSLVLSISAP NLMYLYLDGE FKDIFLENTP KLVAISVSMY MHEDVTD FEQ SSDYNLVKFL GGVPLLEKLV GYIYFTKYLS IGDDPGR LPL TYIHLKTIEL YQVCFEDA VLVLLRLVTH SPNLKELKVS ASPVQPFPLE EEGFDL FERD YFDYKLPSLE SVKITDASGI RYELEFIRFL LGTSPVLETV TVSSSLSDKD AKMDMVI ELL RYPRVSPRAQ LLFLQD
Source	Mammalian cell
Target Names	At1g13570
Protein Names	Recommended name: F-box/FBD/LRR-repeat protein At1g13570
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.