



Recombinant Human Ribosome biogenesis protein WDR12 (WDR12)

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP872413HU |
| Abbreviation | WDR12 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9GZL7 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | AQLQTRFYT DNKKYAVDDV PFSIPAASEI ADLSNIINKL LKDKNEFHKH VEFDFLIKGQ FLRMPLDKHM EMENISSEEV VEIEYVEKYT APQPEQCMFH DDWISSIKGA EEWILTGSYD KTSRIWSLEG KSIMTIVGHT DVVKDVAWVK KDSLSCLLLS ASMDQTILLW EWNVERNKVK ALHCCRGHAG SVDSIAVDGS GTKFCSGSWD KMLKIWSTVP TDEEDEMEES TNRPRKKQKT EQLGLTRTPI VTLSGHMEAV SSVLWSDAEE ICSASWDHTI RVWDVESGSL KSTLTGNKVF NCISYSPLCK RLASGSTDRH IRLWDPRTKD GSLVLSLSTS HTGWVTSVKW SPTHEQQLIS GSLDNIVKLW DTRSCKAPLY DLAAHEDKVL SVDWTDGLL LSGGADNKLY SYRYSPTTSH VGA |
| Source | E.coli |
| Target Names | WDR12 |
| Protein Names | Recommended name: Ribosome biogenesis protein WDR12 Alternative name(s): WD repeat-containing protein 12 |
| Expression Region | 2-423 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.