



Recombinant Human Ribosome biogenesis protein WDR12 (WDR12)

Product Code	CSB-MP872413HU
Abbreviation	WDR12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZL7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AQLQTRFYT DNKKYAVDDV PFSIPAASEI ADLSNIINKL LKDKNEFHKH VEFDFLIKGQ FLRMPLDKHM EMENISSEEV VEIEYVEKYT APQPEQCMFH DDWISSIKGA EEWILTGSYD KTSRIWSLEG KSIMTIVGHT DVVKDVAWVK KDSLSCLLLS ASMDQTILLW EWNVERNKVK ALHCCRGHAG SVDSIAVDGS GTKFCSGSWD KMLKIWSTVP TDEEDEMEES TNRPRKKQKT EQLGLTRTPI VTLSGHMEAV SSVLWSDAEE ICSASWDHTI RVWDVESGSL KSTLTGNKVF NCISYSPLCK RLASGSTDRH IRLWDPRTKD GSLVSLSLTS HTGWVTSVKW SPTHEQQLIS GSLDNIVKLW DTRSCAPLY DLAAHEDKVL SVDWTDGLL LSGGADNKLY SYRYSPTTSH VGA
Source	Mammalian cell
Target Names	WDR12
Protein Names	Recommended name: Ribosome biogenesis protein WDR12 Alternative name(s): WD repeat-containing protein 12
Expression Region	2-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.