



Recombinant Human Neuroblastoma breakpoint family member 3 (NBPF3)

Product Code	CSB-BP872422HU
Abbreviation	NBPF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H094
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPLTPTVQGF QWTLRGPDVE TSPFGAPRAA SHGVGRHQEL RDPTVPGPTS SATNVSMVVS AGPWSGEKAE MNILEINKKS RPQLAENKQQ FRNLKQKCLV TQVAYFLANR QNNYDYEDCK DLIKSMRLRDE RLLTEEKLAE ELGQAEELRQ YKVLVHSQER ELTQLREKLQ EGRDASRSLN QHLQALLTPD EPDNSQGRDL REQLAEGCRL AQHLVQKLSP ENDDDEDEDV KVEEAKEVQE LYAPREVQKA EEKEVPEDSL EECAITCSNS HHPCESNQPY GNTRITFEED QVDSTLIDSS SHDEWLDAVC IIPENESDHE QEEKGPVSP RNLQESEEEE APQESWDEGD WTLSIPDMS ASYQSDRSTF HSVEEQQVGL ALDIGRHWCD QVKKEDQEAT SPRLSRELLD EKEPEVLQDS LDRFYSTPFE YLELPDLCQP YRSDFYSLQE QHLGLALDLD RMKKDQEEEE DQGPPCPRLS RELPEVVEPE DLQDSLDRWY STPFSYPELP DSCQPYGSCF YSLEEEHVGF SLDVDEIEKY QEGEEDQKPP CPRLNEVLME AEEPEVLQDS LDRCYSTTST YFQLHASFQQ YRSAFYSFEE QDVSLALDVD NRFFTLTVIR HHLAFQMGVI FPH
Source	Baculovirus
Target Names	NBPF3
Protein Names	Recommended name: Neuroblastoma breakpoint family member 3 Alternative name(s): Protein AE2 Protein SHIIIa4
Expression Region	1-633
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.