



# Recombinant Human Neuroblastoma breakpoint family member 3 (NBPF3)

<b>Product Code</b>	CSB-EP872422HU-B
<b>Abbreviation</b>	NBPF3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H094
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPLTPTVQGF QWTLRGPDVE TSPFGAPRAA SHGVGRHQEL RDPTVPGPTS SATNVSMVVS AGPWSGEKAE MNILEINKKS RPQLAENKQQ FRNLKQKCLV TQVAYFLANR QNNYDYEDCK DLIKSMRLRDE RLLTEEKLAELGQAEELRQ YKVLVHSQER ELTQLREKLQ EGRDASRSLN QHLQALLTPD EPDNSQGRDL REQLAEGCRL AQHLVQKLSP ENDDDEDEDV KVEEAQKQVE LYAPREVQKA EEKEVPEDSL EECAITCSNS HHPCESNQPY GNTRITFEED QVDSTLIDSS SHDEWLDVAVC IIPENESDHE QEEKGPVSP RNLQSEEEEE APQESWDEGD WTLSIPDMS ASYQSDRSTF HVEEQVGL ALDIGRHWCD QVKKEDQEAT SPRLSRELLD EKEPEVLQDS LDRFYSTPFE YLELPDLCQP YRSDFYSLQE QHLGLALDLD RMKKDQEEEE DQGPPCPRLS RELPEVVEPE DLQDSLDRWY STPFSYPELP DSCQPYGSCF YSLEEEHVGF SLDVDEIEKY QEGEEDQKPP CPRLNEVLME AEEPEVLQDS LDRCYSTTST YFQLHASFQQ YRSYFYSFEE QDVSLALDVD NRFFTLTVIR HHLAFQMGVI FPH
<b>Source</b>	E.coli
<b>Target Names</b>	NBPF3
<b>Protein Names</b>	Recommended name: Neuroblastoma breakpoint family member 3 Alternative name(s): Protein AE2 Protein SHIIIa4
<b>Expression Region</b>	1-633
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.