



# Recombinant Human EH domain-containing protein 4 (EHD4)

<b>Product Code</b>	CSB-EP872437HU-B
<b>Abbreviation</b>	EHD4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H223
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFSWMGRQAG GRERAGGADA VQTVTGGGLRS LYLRKVLPLE EAYRFHEFHS PALEDADFEN KPMILLVGQY STGKTTFIRY LLEQDFPGMR IGPEPTTDSF IAVMYGETEG STPGNALVVD PKKPFRKLSR FGNAFLNRFM CSQLPNQVLK SISVIDSPGI LSGEKQRISR GYDFCQVLQW FAERVDRILL LFDHAKLDIS DEFSEAIKAF RGQDDKIRVV LNKADQVDTQ QLMRVYGALM WSLGKVINTP EVLRVYIGSF WAQPLQNTDN RRLFEAEAQD LFRDIQSLPQ KAAVRKLNLD IKRARLAKVH AYIISYLKKE MPSVFGKENK KRELISRLPE IYIQLQREYQ ISAGDFPEVK AMQEQLENYD FTKFHSLKPK LIEAVDNMLS NKISPLMNLISQEETSTPTQ LVQGGAFDGT TEGPFNQGYG EGAKEGADEE EWVVAKDKPV YDELFTLSP INGKISGVNA KKEMVTSKLP NSVLGKIWKL ADCDCDGM LD EEFALAKHL IKIKLDGYEL PSSLPPHLPV PSHRKSPLKA D
<b>Source</b>	E.coli
<b>Target Names</b>	EHD4
<b>Protein Names</b>	Recommended name: EH domain-containing protein 4 Alternative name(s): Hepatocellular carcinoma-associated protein 10/11 PAST homolog 4
<b>Expression Region</b>	1-541
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.