



Recombinant Human Matrix metalloproteinase-27 (MMP27)

Product Code	CSB-MP872482HU
Abbreviation	MMP27
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H306
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	YG YTLPGWRKYN LTYRIINYTP DMARAAVDEA IQEGLEVWSK VTPLKFTKIS KGIADIMIAF RTRVHGRCPR YFDGPLGLV LG HAFPPGPGLG GDTHFDEDEN WTKDGAGFNL FLVAAHEFGH ALGLSHSNDQ TALMFPNYVS LDPRKYPLSQ DDINGIQSIY GGLPKPAKP KEPTIPHACD PDLTFDAITT FRREVMFFKG RHLWRIYYDI TDVEFELIAS FWPSLPADLQ AAYENPRDKI LVFKDENFWM IRGYAVLPDY PKSIHTLGFP GRVKKIDAAV CDKTRRKYF FVGIWCWRFD EMTQTMCKGF PQRVVKHFPG ISIRVDAAFQ YKGFSSSRG SKQFEYDIKT KNITRIMRTN TWFQCKECPN SSFGFDINKE KAHSGGIKIL YHKSLSLFIF GIVHLLKNTS IYQ
Source	Mammalian cell
Target Names	MMP27
Protein Names	Recommended name: Matrix metalloproteinase-27 Short name= MMP-27 EC= 3.4.24.-
Expression Region	99-513
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP s are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.