



Recombinant Human Peptidase inhibitor R3HDML (R3HDML)

Product Code	CSB-MP872491HU
Abbreviation	R3HDML
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H3Y0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	HISV RDMNALLDYH NHIRASVYPP AANMEYMWVD KRLARAAEAW ATQCIWAHGP SQLMRYVGQN LSIHSGQYRS VVDLMKSWSE EKWHYLF PAP RDCNPHCPWR CDGPTCSHYT QMVWASSNRL GCAIHTCSSI SVWGNTWHRA AYLVCNYAIK GNWIGESPYK MGKPCSSCPP SYQGSCNSNM CFKGLKSNKF TWF
Source	Mammalian cell
Target Names	R3HDML
Protein Names	Recommended name: Peptidase inhibitor R3HDML Alternative name(s): Cysteine-rich secretory protein R3HDML
Expression Region	57-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.