



# Recombinant Human F-box only protein 44 (FBXO44)

<b>Product Code</b>	CSB-EP872497HU-B
<b>Abbreviation</b>	FBXO44
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H4M3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAVGNINELP ENILLELFTH VPARQLLLNC RLVCSLWRDL IDLVTLWKRK CLREGFITED WDQPVADWKI FYFLRSLHRN LLHNPCAEEG FEFWSLDVNG GDEWKVEDLS RDQRKEFPND QVKKYFVTSY YTCLKSQVVD LKAEGYWEEL MDTTRPDIEV KDWFAARPDC GSKYQLCVQL LSSAHAPLGT FQPDPATIQQ KSDAKWREVS HTFSNYPPGV RYIWFQHGGV DTHYWAGWYG PRVTNSSITI GPPLP
<b>Source</b>	E.coli
<b>Target Names</b>	FBXO44
<b>Protein Names</b>	Recommended name: F-box only protein 44 Alternative name(s): F-box protein FBX30 F-box/G-domain protein 3
<b>Expression Region</b>	1-255
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. This protein belongs to the Fbxs class. It is also a member of the NFB42 (neural F Box 42 kDa) family, similar to F-box only protein 2 and F-box only protein 6. Four alternatively spliced transcript variants encoding two distinct isoforms have been found for this gene.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.