



Recombinant Human Uncharacterized protein C20orf132 (C20orf132)

Product Code	CSB-EP872499HU
Abbreviation	C20orf132
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H579
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSSKHRICSQ EEVVIPCAYD SDESVDLEL SNLEIHKKGS SSIELTDLDI PDIPGLHCEP LSHSPRHLTQ QDPLSEAIVE KLIQSIQKVF NGELKGELEK LKFLGDLSSL SQALPYDETA KFSIHSIAD IVHTLNLVQ EERPHSLSSS MRQEVFVTIA DLSYQDVHLL LGSEDRAELF SLTIKSIITL PSVRTLTQIQ EIMPNGTCNT ECLYRQTFQA FSEMLQSLVV KDPHLENLDT IIKLPLRFQR LGHLVALMAL LCGDPQEKVA EEAAEGIHSL LHITLRLKYI THDKKDQQNL KRALTKCREF LELHSSAAKC FYNCPFRIAQ VFEGFLDSNE LCQFIMTTFD TLKTLKHPCI QRSAGELLLT LAKNTESQFE KVPEIMGVIC AQLSIISQPR VRQQIINTVS LFISRPKYTD IVLSFLLCHP VPYNRHLAEV WRMLSVELPS TTWILWRLLR KLQKCHNEPA QEKMAYVAVA VSP
Source	E.coli
Target Names	MROH8
Protein Names	Recommended name: Uncharacterized protein C20orf132
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.