



Recombinant Human PDZ domain-containing protein MAGIX (MAGIX)

Product Code	CSB-EP872508HU-B
Abbreviation	MAGIX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6Y5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEPRTGGAAN PKGSRGSRGP SPLAGPSARQ LLARLDARPL AARAAVDVAA LVRRAGATLR LRRKEAVSVL DSADIEVTDS RLPHATIVDH RPQHRWLETC NAPPQLIQGK AHSAPKPSQA SGHFSVELVR GYAGFGLTLG GGRDVAGDTP LAVRGLLDKG PAQRGCRLEV GDVVLHINGE STQGLTHAQA VERIRAGGPQ LHLVIRPLE THPGKPRGVG EPRKGVVPSW PDRSPDPGGP EVTGSRSSST SLVQHPPSRT TLKKTGRSPE PSPEAAADGP TVSPERRAE DPNDQIPGSP GPWLVPSEER LSRALGVRGA AQFAQEMAAG RRRH
Source	E.coli
Target Names	MAGIX
Protein Names	Recommended name: PDZ domain-containing protein MAGIX
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.