



Recombinant Human Transducin-like enhancer protein 6 (TLE6)

Product Code	CSB-EP872514HU-B
Abbreviation	TLE6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H808
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MTSRDQPRPK GPPKSTSPCP GISNNESSPT LNYQGILNRL KQFPRFSPHF AAELESIYYS LHKIQQDVAE HHKQIGNVLQ IVESCSQLQG FQSEEVSPA PASPGTPQQV KDKTLQESSF EDIMATRSSD WLRRLGEDN QPETQLFWDK EPFWHDTLT EQLWRIFAGV HDEKAKPRDR QQAPGLGQES KAPGSCDPGT DPCPEDASTP RPPEASSSPP EGSQDRNTSW GVVQEPPGRA SRFLQISWD PEDFEDAWKR PDALPGQSKR LAVPCKLEKM RILAHGELVL ATAISSFTRH VFTCGRRGIK VWSLTGQVAE DRFPESHLP QTPGAFLRTC LLSSNSRSL TGGYNLASVS VWDLAAPSLH VKEQLPCAGL NCQALDANLD ANLAFASFTS GVVRIWDLRD QSVVRDLKGY PDGVKSIVVK GYNIWTGGPD ACLRCWDQRT IMKPLEYQFK SQIMSLSHSP QEDWVLLGMA NGQQWLQSTS GSQRHVMVGQK DSVILSVKFS PFGQWWASVG MDDFLGVYSM PAGTKVFEVP EMSPVTCCDV SSNNRLVVTG SGEHASVYQI TY </p>
Source	E.coli
Target Names	TLE6
Protein Names	Recommended name: Transducin-like enhancer protein 6
Expression Region	1-572
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.