



Recombinant Human Actin-related protein 8 (ACTR8)

Product Code	CSB-BP872523HU
Abbreviation	ACTR8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H981
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTQAEKGDTE NGKEKGGEKE KEQRGV KRPI VPALVPESLQ EQIQSNFIIV IHPGSTLLRI GRATDTLPAS IPHVIARRHK QGQPLYKDS WLLREGLNKP ESNEQRQNGL KMVDQAIWSK KMSNGTRRIP VSPEQARSYN KQMRPAILDH CSGNKWTNTS HHPEYLVGEE ALYVNPLDCY NIHWPIRRGQ LNIHPGPGGS LTAVLADIEV IWSHAIQKYL EIPLKDLKYY RCILLIPDIY NKQHVKELVN MILMKMGFSG IVVHQESVCA TYGSGLSSTC IVDVGDQKTS VCCVEDGVSH RNTRLCLAYG GSDVSRCFYW LMRAGFPYR ECQLTNKMDC LLLQHLKETF CHLDQDISGL QDHEFQIRHP DSPALLYQFR LGDEKLQAPM ALFYPTFGI VGQKMTTLQH RSQGDPEDPH DEHYLLATQS KQE QSAKATA DRKSASKPIG FEGDLRGQSS DLPERLHSQE VDLGSAQGDG LMAGNDSEEA LTALMSRKTA ISLFEGKALG LDKAILHSID CCSSDDTKKK MYSSILVGG GLMFHKAQEF LQHRILNKMP PSFRRIENV DVITRPKDM D PRLIAWKGGA VLACLDTTQE LWIYQREWQR FGVRMLRERA AFVW
Source	Baculovirus
Target Names	ACTR8
Protein Names	Recommended name: Actin-related protein 8 Short name= hArp8 Alternative name(s): INO80 complex subunit N
Expression Region	1-624
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.