



Recombinant Human Ankyrin repeat family A protein 2 (ANKRA2)

Product Code	CSB-EP872524HU
Abbreviation	ANKRA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H9E1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDTSTNLDIG AQLLIVEECPS TYSLTGMPDI KIEHPLDPNS EEGSAQGVAM GMKFILPNRF DMNVCSRFBK SLNEEDSKNI QDQVNSDLEV ASVLFKAECN IHTSPSPGIQ VRHVYTPSTT KHFSPIKQST TLTNKHARGNE VSTTPLLANS LSVHQLAAQG EMLYLATRIE QENVINHTDE EGFTPLMWAA AHGQIAVVEF LLQNGADPQL LGKGRESALS LACSKGYTDI VKMLLDGCGVD VNEYDWNGGT PLLYAVHGNH VKCVKMLLES GADPTIETDS GYNSMDLAVA LGYRSVQQVI ESHLLKLLQN IKE
Source	E.coli
Target Names	ANKRA2
Protein Names	Recommended name: Ankyrin repeat family A protein 2 Alternative name(s): RFXANK-like protein 2
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.