



Recombinant Human COP9 signalosome complex subunit 7b (COPS7B)

Product Code	CSB-EP872526HU-B
Abbreviation	COPS7B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H9Q2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AGEQKPSSN LLEQFILLAK GTSGSALTAL ISQVLEAPGV YVFGELLELA NVQELAEGAN AAYLQLLNLF AYGTYPDYIA NKESSLPELST AQQNKLKHLT IVSLASRMKC IPYSVLLKDL EMRNLRELED LIIEAVYTDI IQGKLDQRNQ LLEVDFCIGR DIRKKDINNI VKTLHEWCDG CEAVLLGIEQ QVLRANQYKE NHNRTQQQVE AEVTNIKKTL KATASSSAQE MEQQLAEREC PPHAEQRQPT KKMSKVKGLV SSRH
Source	E.coli
Target Names	COPS7B
Protein Names	Recommended name: COP9 signalosome complex subunit 7b Short name= SGN7b Short name= Signalosome subunit 7b Alternative name(s): JAB1- containing signalosome subunit 7b
Expression Region	2-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.