



Recombinant Human Chromatin modification-related protein MEAF6 (MEAF6)

Product Code	CSB-MP872530HU
Abbreviation	MEAF6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HAF1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAMHNKAAPP QIPDTRRELA ELVKRKQELA ETLANLERQI YAFEGSYLED TQMYGNIIRG WDRYLTNQKN SNSKNDRNR KFKEAERLFS KSSVTSAAAV SALAGVQDQL IEKREPGSGT ESDTSPDFHN QENEPSQEDP EDLDGSGVQGV KPQKAASSTS SGSHHSSHKK RKNKNRHRID LKLNKKPRAD Y
Source	Mammalian cell
Target Names	MEAF6
Protein Names	Recommended name: Chromatin modification-related protein MEAF6 Short name= MYST/Esa1-associated factor 6 Alternative name(s): Esa1-associated factor 6 homolog Short name= Protein EAF6 homolog Short name= hEAF6 Sarcoma anti
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.