



Recombinant Human Homeobox protein SEBOX (SEBOX)

Product Code	CSB-EP872535HU-B
Abbreviation	SEBOX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HB31
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGGSGVGTAW HGPLARPSGT LPFASSMPSP VDASSADGGS GLGSHRRKRT TFSKGQLEL ERAFAAWPYP NISTHEHLAW VTCLPEAKVQ VWFQKRWAKI IKNRKSGILS PGSECPQSSC SLPDTLQQPW DPQMGPQPPP SSGTPQRTSV CRHSSCPAPG LSPRQGWEGA KAVAPWGSAG ASEVHPSLER ATPQTSLGSL SDLIYALAI VNV DHS
Source	E.coli
Target Names	SEBOX
Protein Names	Recommended name: Homeobox protein SEBOX Alternative name(s): Homeobox OG-9 Skin-, embryo-, brain- and oocyte-specific homeobox
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.