



Recombinant Human Band 4.1-like protein 5 (EPB41L5)

Product Code	CSB-YP872549HU
Abbreviation	EPB41L5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HCM4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> MLSFFRRTLK RRSRMRKHAEK ERLREAQRAA THIPAAGDSK SIITCRVSLI DGTDVSVLDP KKAKGQELFD QIMYHLDLIE SDYFGLRFMD SAQVAHWLDG TKSIIKQVKI GSPYCLHLRV KFYSSEPNNL REELTRYLFV LQLKQDILSG KLDCPFDTAV QLAAYNLQAE LGDYDLAEHS PELVSEFRFV PIQTEEMELA IFEKWKEYRG QTPAQAEYNY LNKAKWLEMY GVDMHVVKAR DGNDYSLGLT PTGVLVFEGD TKIGLFFWPK ITRLDFKKNK LTLVVVEDDD QGKEQEHTFV FRLDHPKACK HLWKCAVEHH AFFRLRGPVQ KSSHRSGFIR LGSFRFRYSK TEYQTTKTKN ARRSTSFERR PSKRYSRRTL QMKACATKPE ELSVHNNVST QSNGSQQAWG MRSALPVSPS ISSAPVPVEI ENLPQSPGTD QHDKRCIPLN IDLLNSPDLL EATIGDVIGA SDTMETSQAL NDVNVATRLP GLGEPEVEYE TLKDTSEKLN QLEMENSPLL SPRSNIDVNI NSQEEVVKLT EKCLNNVIES PGLNVMRVPP DFKSNILKAQ VEAHVHVTKE DSLLSHKNAN VQDAATNSAV LNENNVPLPK ESLETLMLIT PADSGSVLKE ATDELDAALL SLTENLIDHT VAPQVSSTSM ITPRWIVPQS GAMSNGLAGC EMLLTGKEGH GNKDGISLIS PPAPFLVDAV TSSGPILAE AVLKQKCLLT TEL </p>
Source	Yeast
Target Names	EPB41L5
Protein Names	Recommended name: Band 4.1-like protein 5
Expression Region	1-733
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.