



# Recombinant Halobacterium salinarum Signal recognition 54 kDa protein (srp54)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP872625HTM-B  |
| <b>Abbreviation</b>      | srp54  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q9HMN5   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Halobacterium salinarum (strain ATCC 700922 / JCM 11081 / NRC-1)<br>(Halobacterium halobium)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | MVLDDLGSSL RGTLDTLRGK SRISEEDVDK VVKEIQRSL QADVVDLVM<br>DLSESITERS LEEPPGGTS ARDHVLRIVY EELVGLVGDS TPIPLEEQTI<br>MLAGLQGS GK TTTAAKMAWW FSKKGLRPAI VQTDTRPGA YEQAEMADR<br>AEVDFYGE PD SDDPVQIARD GLAATADADV HIVDTAGRHA LEDDLIAEIE<br>EIEGVVEPDR SLLVLDAAIG QGAKQQARQF DASIGIDGVA ITKLDGTAKG<br>GGSLTAVNET DSTIAFLGSG ETVQDIERFE PDSFISRL LG MGD LKQLTER<br>VERAMEETQD DDEDWDPE DL MKGEFTLHDM RKQMAMDNM<br>GPLDQVMDMI PGMGGGLMDE LPDDAMDVTQ ERMRGFEVIM<br>DSMTEAMAD PRSVGASQIR RIARGSGQPE DRVRELLDQH<br>KMMAQTMQQF QGMGDGDMQR MMSKMQQQGG GGGGGFGGMF |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | srp54  |
| <b>Protein Names</b>     | Recommended name: Signal recognition 54 kDa protein Alternative name(s):<br>SRP54  |
| <b>Expression Region</b> | 1-460  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.   |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.