



# Recombinant Mouse Serine/threonine-protein kinase 19 (Stk19)

<b>Product Code</b>	CSB-MP872828MO
<b>Abbreviation</b>	Stk19
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9JHN8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNRKRLLAS EAFGVKRRRA PGPVRADPLR TRAGSAREAI EELVKLFPRG LFEDALPPIA LRSQVYSLVP DRTVADLQLK ELQELGEIRI IQLGFDLDAH GIVFTEDYRT RVLKACDGRP CAGAVQKFLA SVLPACGDLS FQQDQMTQTY GFRDPEITQL VNAGVLTVRD AGSWWLAVPG AGRFIKCFVK GRQAVLSMVR KAKYRELALS ELLGRRAPLA VRLGLAYHVH DLIGAQLVDC VPTTSGTLLR LPDT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Stk19
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase 19 EC= 2.7.11.1 Alternative name(s): Protein RP1
<b>Expression Region</b>	1-254
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a serine/threonine kinase which localizes predominantly to the nucleus. Its specific function is unknown; it is possible that phosphorylation of this protein is involved in transcriptional regulation. This gene localizes to the major histocompatibility complex (MHC) class III region on chromosome 6 and expresses two transcript variants.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.