



Recombinant Rat Cullin-5 (Cul5)

Product Code	CSB-YP872849RA
Abbreviation	Cul5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9JJ31
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	ATSNLLKNK GSLQFEDKWD FMRPIVLKLL RQESVTKQQW FDLFSDVHAV CLWDDKGSSK IHQALKEDIL EFIKQAQARV LSHQDDTALL KAYIVEWRKF FTQCDILPKP FCQLEVTLG KQSSNKSNM EDSIVRKLML DTWNESIFSN IKNRLQDSAM KLVHAERLGE AFDSQLVIGV RESYVNLCSN PEDKLQIYRD NFEKAYLDST ERFYRTQAPS YLQQNGVQNY MKYADAKLKE EEKRALRYLE TRRECNSVEA LMECCVNALV TSFKETILAE CQGMIKRNET EKLHLMFSLM DKVPGGIEPM LKDLEEHIS AGLADMVAAA ETITTDSEKY VEQLLTLFNR FSKLVKEAFQ DDPRLTARD KAYKAVVND TIFKLELPLK QKGVGLKTQP ESKPELLAN YCDMLLRKTP LSKKLTSEEI EAKLKEVLLV LKYVQNKDVF MRYHKAHLTR RLILDISADS EIEENMVEWL REVGMPADYV NKLARMFQDI KVSEDLNQAF KEMHKNNKLA LPADSVNIKI LNAGAWSRSS EKVFVSLPTE LEDLIPEVEE FYKKNHSGRK LHWHLMSNG IITFKNEVGQ YDLEVTTFQL AVLFAWNQRP REKISFENLK LATELPDAEL RRTLWLSLVA FPKLKRQVLLY DPQVNSPKDF TEGTLFSVNQ DFSLIKNAKV QKRGKINLIG RLQLTTERMR EEENEGIVQL RILRTQEAI QIMKMRKKIS NAQLQTELVE ILKNMFLPQK KMIKEQIEWL IEHKYIRRDE ADINTFIYMA
Source	Yeast
Target Names	Cul5
Protein Names	Recommended name: Cullin-5 Short name= CUL-5 Alternative name(s): Vasopressin-activated calcium-mobilizing receptor 1 Short name= VACM-1
Expression Region	2-780
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.