



# Recombinant Rat Apoptosis-inducing factor 1, mitochondrial (Aifm1)

<b>Product Code</b>	CSB-YP872890RA
<b>Abbreviation</b>	Aifm1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9JM53
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GLSPEEKQR RAIASAAEGG SVPPIRVPSH VPFLDIGGGT AAFAAARSIR ARDPGARVLI VSEDELPYM RPPLSKELWF SDDPNVTKTL QFRQWNGKER SIYFQPPSFY VSAQDLPHIE NGGVAULTGK KVVHLDVRGN MVKLNDSQI TFEKCLIATG GTPRSLSAID RAGAEVKSRT TLFKIGDFR ALEKISREVK SITVIGGGFL GSELACALGR KSQASGIEVI QLFPEKGNMG KILPEYLSNW TMEKVKREGV KVMPNAIVQS VGVSGGKLLI KLDGRKVVET DHIVTAVGLE PNVELAKTGG LEIDSDFGGF RVNAELQARS NIWVAGDAAC FYDIKLGRRR VEHHDHAVVS GRLAGENMTG AAKPYWHQSM FWSDLGPDVG YEAIGLVDSS LPTVGVFAKA TAQDNPKSAT EQSGTGIRSE SETESEASEI TIPSPDPAVP QVPVEGEDYG KGVIFYLRDK VVGIVLWNV FNRMPIARKI IKDGEQHEDL NEVAKLFNIH ED
<b>Source</b>	Yeast
<b>Target Names</b>	Aifm1
<b>Protein Names</b>	Recommended name: Apoptosis-inducing factor 1, mitochondrial EC= 1.-.-.- Alternative name(s): Programmed cell death protein 8
<b>Expression Region</b>	102-612
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a flavoprotein essential for nuclear disassembly in apoptotic cells that is found in the mitochondrial intermembrane space in healthy cells. Induction of apoptosis results in the translocation of this protein to the nucleus where it effects chromosome condensation and fragmentation. In addition, this gene product induces mitochondria to release the apoptogenic proteins cytochrome c and caspase-9. Several alternative transcripts encoding different isoforms have been identified for this gene.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.