



# Recombinant Mouse Sialidase-3 (Neu3)

<b>Product Code</b>	CSB-YP872896MO
<b>Abbreviation</b>	Neu3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9JMH7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEEVPPYSLS STLFQQEEQS GVTYRIPALL YLPPTHFLA FAEKRTSVRD EDAACLVLRG GLMKGRSVQW GPQRLLMEAT LPGHRTMNPC PVWEKNTGRV YLFFICVRGH VTERCQIVWG KNAARLCFLC SEDAGCSWGE VKDLTEEVIG SEVKRWATFA VGPGHGIQLH SGRLIIPAYA YYVSRWFLCF ACSVKPHSLM IYSDDFGVTW HHGKFIEPQV TGECQVAEVA GTAGNPVLYC SARTPSRFRA EAFSTDSGGC FQKPTLNPQL HEPRTGCQGS VVSFRPLKMP NTYQDSIGKG APATQKCP LL DSPLEVEKGA ETPSATWLLY SHPTSKRKRI NLGIYYNRNP LEVNCWSRPW ILNRGPSGYS DLAVVEEQDL VACLFECEGK NEYERIDFCL FSDHEVLSCE DCTSPSSD
<b>Source</b>	Yeast
<b>Target Names</b>	Neu3
<b>Protein Names</b>	Recommended name: Sialidase-3 EC= 3.2.1.18 Alternative name(s): Ganglioside sialidase Membrane sialidase N-acetyl-alpha-neuraminidase 3
<b>Expression Region</b>	1-418
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene product belongs to a family of glycohydrolytic enzymes which remove sialic acid residues from glycoproteins and glycolipids. It is localized in the plasma membrane, and its activity is specific for gangliosides. It may play a role in modulating the ganglioside content of the lipid bilayer.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.