



Recombinant Bacillus halodurans Phosphatidylserine decarboxylase proenzyme (psd)

| | |
|--------------------------|---|
| Product Code | CSB-YP873066BQS |
| Abbreviation | psd |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9KDA3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bacillus halodurans (strain ATCC BAA-125 / DSM 18197 / FERM 7344 / JCM 9153 / C-125) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MRKTFFRLCF ELSSHPVLSS MLQAFTTСКА SRWLIPSFVR VYNINGQAE KPLHTYQSLQ EVFTRTLTEN CRPIDLSPKS IVSPVDGVA EQGTLSEAN FVVKNQTYTL EEMLGKKEKA KLYREGTYLL FYLSPSHYHR IHSPVNGTIK EQWTLGNKSA PVNNLGLRYG KRPLSRNYRL LTELEAEEGR CIVAKIGALN VNSIVPTHQS EHVDPGEEIG YFAFG |
| Source | Yeast |
| Target Names | psd |
| Protein Names | Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain |
| Expression Region | 1-225 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.