



Recombinant Arabidopsis thaliana Putative F-box protein At3g24700 (At3g24700)

Product Code	CSB-MP873278DOA
Abbreviation	At3g24700
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LJ39
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MLTDLPLDLE SEILSRVPAT SLQRLKTTCK RWYALFRDPR FVKKNLGKAA THVIFDNRSG YSMTDINSLI HSINLRGIQN SFDPSIGVDV KLNVLKDPRH DKISHIISHC DGLLLCKTED YGRLVVWNPC TGQIKWIQAN NMLMDVYVLG YVNNNKSCNS YKILNFGILP LNSSHDKNSK IYEFNSDSWR ILDHVSPGYF AISKAMTLKG NAYWFASDWS GTKTRTRIYF AKEGIQVHQE IAQKPKKCGP FLVSYVPSLV QIQSGNKQSS
Source	Mammalian cell
Target Names	At3g24700
Protein Names	Recommended name: Putative F-box protein At3g24700
Expression Region	1-270
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.