



# Recombinant Arabidopsis thaliana Putative protease Do-like 12, mitochondrial (DEGP12)

<b>Product Code</b>	CSB-BP873288DOA
<b>Abbreviation</b>	DEGP12
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LK70
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SRIATI VLPFALTRGR KIHTMSKDEE WWKKIRKSPV VDELMLESVV EVFTDSTKYS KVKPWQTLNQ ESYGGSGFAI AGKKILTNAH VVEGMNDHIF VHVKRHGSQV KYKAKVQKIA HECDLAILEI DSDEFWKGMM PLEFGDIPPL NEIVYVVGYP KAGETICVTK GVVTVGKTGN YLRSSTKLLT IHIDATTYGG NSGGPVITGD KVLGVLFQIL GDKKSTGVVI PTPIIRHFIT GAEESSHNAV FGSLVLSCQS MKNAQIRNHF KMSPETTGIL INKINSSSGA HKILRKDDII LAIDGVPVLS EMRRISFNHF ISMKKPDENI LVKVLKRGKE HEYNISLKPV KPHIQVQQYY NLPSYYIFGG FVFPVPLTKSY IDDKYYKITD EQHVVISQVM PDDINKGYSN FKDLQVEKVN GVKVKNLHL RELIEGCFSSK DLRLDLENDK VMVLNYESAK KATFEILERH NIKSAWASE
<b>Source</b>	Baculovirus
<b>Target Names</b>	DEGP12
<b>Protein Names</b>	Recommended name: Putative protease Do-like 12, mitochondrial EC= 3.4.21.-
<b>Expression Region</b>	25-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.