



Recombinant Arabidopsis thaliana Probable carboxylesterase 1 (CXE1)

Product Code	CSB-MP873296DOA
Abbreviation	CXE1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LMA7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDSEIAFDYS PRFRIFKNGG IERLVPETFV PPSLNPENGV VSKDAVYSPE KNLSLRILYP QNSVYETGEK KIPLLVFYFHG GGFIMETAFS PIYHTFLTSA VSATDCIAVS VEYRRAPEHP IPTLYEDSWD AIQWIFTHIT RSGPEDWLNK HADFSKVFLA GDSAGANIAH HMAIRVDKEK LPPENFKISG MILFHPYFLS KALIEEM EVE AMRYERLWR IASPDGNGV EDPWINVVG S DLTGLGCRRV LVMVAGNDVL ARGGWSYVAE LEKSGWIGKV KVMETKEEGH VFHLRDPDSE NARRVLRNFA EFLKEETF
Source	Mammalian cell
Target Names	CXE1
Protein Names	Recommended name: Probable carboxylesterase 1 Alternative name(s): AtCXE1 EC= 3.1.1.1
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.