



Recombinant Arabidopsis thaliana F-box protein At1g59680 (At1g59680)

Product Code	CSB-EP873313DOA-B
Abbreviation	At1g59680
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LQ46
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTTMSDLSVD LVGEILSRVP LTSLSAVRCT CKSWNTLSKH QIFGKAELAA TKQFLGFTVM DYKVCSLRFD LQGIRNDGDD FVDHGEVYDV FHSDGLLLCV TKDHWRLVW NPYLGQIRWI RARKQFGYSN MFSLGYNNSN RNHKIIEVSY VYDKTACPQY SEIYDFNSSS WRLVEINPVW FLHSDRVSLK GNTYFFAQDL LKDAKIENYL LCFDFTAERF GPHLPLPVHS YDDVTLSQVVR QEQLALLYGN NETDDSLEIW ITNQIDPNAV SWSIFLKVDI KPLIGFPKDF DPGSFFIDEE KKVVVVYDLD GHLYSTKKNC AYQRLYIIGE DGYFTSVKIK CSDSPEFSAY APSLVITYKPT NYAK
Source	E.coli
Target Names	At1g59680
Protein Names	Recommended name: F-box protein At1g59680
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.