



Recombinant Arabidopsis thaliana Putative F-box protein At3g16590 (At3g16590)

Product Code	CSB-MP873365DOA
Abbreviation	At3g16590
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LUS5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MPTKLPLELE DEILLRVPL SLTRFRTVCK RWNTLFNDQR FINNHLACVR PQFILRTEKD SKIYSIGINI DDSLEVRELN LETQGPNKKL KVYRNLFYCD GFLLCALLD EVAVWNPWLR KQTKWIEPKR SRFNLYGLGY DNRRPEKCYK ILGFGYGYSS EINGSYNRIN PRVSVFET NAWKDLKFGF FDWHLRSPRT VLSLNGTLYW IAVRCESGGD GFIQSFDFSR EMFEPFCLLP CKNDFGDTQI LEVFRGDRLS VLEQCPTTNK IKIWWTKNKI SGDRKELVSW RLLMTVSIPN FPRLQDLYSN SQPSYFMDNN DDKRLIVCTC DESGKPCIYI VKGDRFKKIQ MGFEVEPWPF HLVYVPSLVP IPLA
Source	Mammalian cell
Target Names	At3g16590
Protein Names	Recommended name: Putative F-box protein At3g16590
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.