



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At3g17280 (At3g17280)

Product Code	CSB-MP873366DOA
Abbreviation	At3g17280
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LUU3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTTISDLPYD LLPEILSRP TKSIPKLKTT CKKWYALFKD PKFVEKKLGK AARETVFLMN HEVNSISVDI HGIPKGYSVS MDFTGTLTIP EGS DLEIFRI HHCNGLFLCA TMNCRLVWN PCTGQITWII PRTRYDSDDI YALGCGDDKS SSLHSYKILR CCDDNQQKPV SEIYDFSSSS WRVLDGVTAN CFIECNGVAL KESAYWYASD KRETPKGKFI LRFDFATERF ARLCLPLNFQ RDRDNKSVVV SVVGEEKLAL LQQFDHRVHG LKYSKIKIWW TDTKIGEGKD LSWSNILVEE LADDNLPSVT SFLLDEEKKV AVCSDAVCSD TDEDED RRR IYIVGEGVDE FVYDEVSTET SHNWPFLVSY VPNLVHIEKD APIVEV
Source	Mammalian cell
Target Names	At3g17280
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At3g17280
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.