



Recombinant Arabidopsis thaliana PsbP domain-containing protein 6, chloroplastic (PPD6)

Product Code	CSB-MP873397DOA
Abbreviation	PPD6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LXX5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	REV EVGSYLPLSP SDPSFVLFKA KPSDTPALRA GNVQPYQFVL PPNWKQLRIA NILSGNYCQP KCAEPWIEVK FENEKQGKVQ VVASPLIRLT NKPNIATIEDL GEPEKVIASL GPFVTGNSYD SDELLKTSIE KIGDQTYKY VLETPFALTG SHNLAKATAK GSTVVLFVVS ATEKQWQSSQ KTLEAILDSF QL
Source	Mammalian cell
Target Names	PPD6
Protein Names	Recommended name: PsbP domain-containing protein 6, chloroplastic Alternative name(s): OEC23-like protein 1
Expression Region	68-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.