



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At5g02980 (At5g02980)

Product Code	CSB-BP873409DOA
Abbreviation	At5g02980
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LYY7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNESSESPPP SRTFSSLPYD IILNCLARVS RYHYPTLSLV SKEFQSLIAS RELYATRSRI GKTERFLYIC LNLTKSNPKY RWFTLPPVPN EQKLLPVPLF TYHLNSSTVS STDSEIYIIG GLVWGNRSKK ASIFDCRSHQ TRRLPKMRFP RASAAAHVID GKIYVIGGGE IRGEVYDPTT QTWLTPVDH TTEECQKVYD KHGVNICFVE IDNLLCQTFV FNGKLYWRHP RGDDFGWARV KGVEQELFRN HLSYVDKSGG GRRVTVWWKS VVVFQCQGLG YSTEFDTEIW CAEISFERRG LKELWGFVEW SRKVFTCKGC DSPTDFVMDS VLVTY
Source	Baculovirus
Target Names	At5g02980
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At5g02980
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.