



Recombinant Arabidopsis thaliana Putative F-box protein At4g17200 (At4g17200)

Product Code	CSB-MP873427DOA
Abbreviation	At4g17200
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M0M7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTTMSDLSPD LVGEILTRVP MTSLISVRCT CKMWNALSKE GIFFKAARKQ FMGFTMMDSR VCSMKFDLQG IRNNEHDFVD PCIKQIAKLD QIEVSKVLQC DGLLLCVGKD NSRLLVWNPY LGQTRFIKPR KRFNKLEWYA LGYDNNHNYK ILMNYDTGHL FGYEIYDFSS DSWRVLDVSP DCHIPFHQQS VSLNGNTYFL AQEKIIVEGE EEEVVDEIEK FLLCFDFTTE RFGPRLPLPF HSDVLETVAI SCVRDDQLAV LYQRFETWEI WVTTKIDPTA VSWSMFLSVD MEPPTGCQFD DEAGSFFIDE ENKVAVFFNS DLFNPRVENR NRRYPTAYII GQDGYFKSAT LREAPDLVKP DPYVYCCPLV WSSYAPSLVQ L
Source	Mammalian cell
Target Names	At4g17200
Protein Names	Recommended name: Putative F-box protein At4g17200
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.