



Recombinant Arabidopsis thaliana SKP1-like protein 13 (ASK13)

Product Code	CSB-BP873439DOA
Abbreviation	ASK13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M1X5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSKMVMLLSS DGESFQVEEA VAVQSQTIAH MIEDDCVANG VPIANVTGVI LSKVIEYCKK HVVSDSPTEE SKDELKKWDA EFMKALEQSS TLFDVMLAAN YLNKDLLDL GCQTVADMIT GKKPDEIRAL LGIENDFTPE EEEEIRKENQ WAFE
Source	Baculovirus
Target Names	ASK13
Protein Names	Recommended name: SKP1-like protein 13 Short name= AtSK13
Expression Region	1-154
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.