



Recombinant Drosophila melanogaster NIF3-like protein 1 (anon-35F/36A)

Product Code	CSB-MP873601DLU
Abbreviation	anon-35F/36A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NK57
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLRRMNTFSG QKLAADVVEL ENFAPTSWAE KWDNVGLLIE PHREKQIKKI LLTNDLTEPV VKEALEKEAE LIISYHPIIF KPLTRITQSH WKERVVAACL ANDIALYSPH TAWDKKSGGV NDWLSKAVNI ISIRPLEPEL GAPPGTGSGR YIETKMELSQ VVESLQKRIR NSVHVALAVG HTPKTLIQSV GICAGSGASL LKGIQADLII TGEMSHHEVL EFTHNNTTVL LCNHSNSERG FLHEFCPILA KSLNEECLVF VSEVDKPLV TVASDINKEL SAFVDVYKST SK
Source	Mammalian cell
Target Names	anon-35F/36A
Protein Names	Recommended name: NIF3-like protein 1 Alternative name(s): Protein anon-35F/36A
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.