



Recombinant Human Mid1-interacting protein 1 (MID1IP1)

Product Code	CSB-EP873607HU
Abbreviation	MID1IP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NPA3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMQICDTYNQ KHSLFNAMNR FIGAVNNMDQ TVMVPSLLRD VPLADPGLDN DVGVEVGGSG GCLEERTPPV PDSGSANGSF FAPSRDMYSH YVLLKSIRND IEWGVLHQPP PPAGSEEGSA WSKSKDILVDL GHLEGADAGE EDLEQQFHYYH LRGLHTVLSK LTRKANILTN RYKQEIGFGN WGH
Source	E.coli
Target Names	MID1IP1
Protein Names	Recommended name: Mid1-interacting protein 1 Alternative name(s): Gastrulation-specific G12-like protein Mid1-interacting G12-like protein Protein STRAIT11499 Spot 14-related protein Short name= S14R Short name= Spot 14-R
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.