



Recombinant Human Calcium-binding protein 2 (CABP2)

Product Code	CSB-EP873608HU
Abbreviation	CABP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NPB3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GNCAKRPWR RGPKDPLQWL GSPPRGSCPS PSSSPKEQGD PAPGVQGYSV LNSLVGPACI FLRPSIAATQ LDRELRP EEI EELQVAFQEF DRDRDGYIGC RELGACMRTL GYMPTE MELI EISQQISGGK VDFEDFVELM GPKLLAETAD MIGVREL RDA FREFDTNGDG RISVGELRAA LKALLGERLS QREVDEILQD VDLNGDGLVD FEEFVRMMSR
Source	E.coli
Target Names	CABP2
Protein Names	Recommended name: Calcium-binding protein 2 Short name= CaBP2
Expression Region	2-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The product of this gene belongs to a subfamily of calcium binding proteins, which share similarity to calmodulin. Calcium binding proteins are an important component of calcium mediated cellular signal transduction. Expression of this gene is retina-specific. The specific function of this gene is unknown. Alternatively spliced transcript variants encoding different isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.