



Recombinant Human Fanconi anemia group F protein (FANCF)

Product Code	CSB-EP873614HU-B
Abbreviation	FANCF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NPI8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MESLLQHLDR FSELLAVSST TYVSTWDPAT VRRALQWARY LRHIHRRFGR HGPIRTALER RLHNQWRQEG GFGRGPVPGI ANFQALGHCD VLLSLRLLN RALGDAARYH LVQQLFPGPG VRDADEETLQ ESLARLARRR SAVHMLRFNG YREPNLQED SLMKTQAELL LERLQEVGKA EAERPARFLS SLWERLPQNN FLKVIALL QPPLSRRPQE ELEPGIHKSP GEGSQVLVHW LLGNSEVFAA FCRALPAGLL TLVTSRHPAL SPVYLGLLTD WGQRLHYDLQ KGIWVGTESQ DVPWEELHNR FQSLCQAPP LKDKVLTAL TCKAQDGD FE VPGLSIWTDL LLALRSGAFR KRQVLGLSAG LSSV
Source	E.coli
Target Names	FANCF
Protein Names	Recommended name: Fanconi anemia group F protein Short name= Protein FACF
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The Fanconi anemia complementation group (FANC) currently includes FANCA, FANCB, FANCC, FANCD1 (also called BRCA2), FANCD2, FANCE, FANCF, FANCG, FANCI, FANCI (also called BRIP1), FANCL, FANCM and FANCN (also called PALB2). The previously defined group FANCH is the same as FANCA. Fanconi anemia is a genetically heterogeneous recessive disorder characterized by cytogenetic instability, hypersensitivity to DNA crosslinking agents, increased chromosomal breakage, and defective DNA repair. The members of the Fanconi anemia complementation group do not share sequence similarity; they are related by their assembly into a common nuclear protein complex. This gene encodes the protein for complementation group F.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.