



Recombinant Human Acetyl-coenzyme A synthetase, cytoplasmic (ACSS2)

Product Code	CSB-EP873629HU-B
Abbreviation	ACSS2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NR19
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGLPEERVRS GSGSRGQEEA GAGGRARSWS PPPEVSRSAH VPSLQRYREL HRRSVEEPRE FWG DIAKEFY WKTPCPGPFL RYNFDVTKGK IFIEWMKGAT TNICYNVLDR NVHEKKLGDK VAFYWEGNEP GETTQITYHQ LLVQVCQFSN VLRKQGIQKG DRVAIYMPMI PELVVAMLAC ARIGALHSIV FAGFSSESLC ERILDSSCSL LITDAFYRG EKLVLNKELA DEALQKCQEK GFPVRCCIVV KHLGRAELGM GDSTSQSPPI KRSCP DVQIS WNQGIDLWWH ELMQEAGDEC EPEWCDAEDP LFILYTS GST GPKPGVVHTV GGYMLYVATT FKYVDFDFAE DVFWCTADIG WITGHSYVTY GPLANGATSV LFEGIPTYPD VNRLWSIVDK YKVTKFYTAP TAIRLLMKFG DEPVTKHSRA SLQVLGTVGE PINPEAWLWY HRVVG AQRPC IVDTFWQTET GGHMLTPLPG ATPMKPGSAT FPFFGVAPAI LNESGEELEG EAEGYLVFKQ PWP GIMRTVY GNHERFETTY FKKFP GYYVT GDGCQRDQDG YYWITGRIDD MLNVSGHLLS TAEVESALVE HEAVAEAAVV GHPHPVKGEC LYCFVTLCDG HTFSPKLTEE LKKQIREKIG PIATPDYIQN APGLPKTRSG KIMRRVLRKI AQNDHDLGDM STVADPSVIS HLFSHRCLTI Q
Source	E.coli
Target Names	ACSS2
Protein Names	Recommended name: Acetyl-coenzyme A synthetase, cytoplasmic EC= 6.2.1.1 Alternative name(s): Acetate--CoA ligase Acetyl-CoA synthetase Short name= ACS Short name= AceCS Acyl-CoA synthetase short-chain family member 2
Expression Region	1-701
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a cytosolic enzyme that catalyzes the activation of acetate for use in lipid synthesis and energy generation. The protein acts as a monomer



and produces acetyl-CoA from acetate in a reaction that requires ATP. Expression of this gene is regulated by sterol regulatory element-binding proteins, transcription factors that activate genes required for the synthesis of cholesterol and unsaturated fatty acids. Alternative splicing results in multiple transcript variants.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.