



Recombinant Human SERTA domain-containing protein 4 (SERTAD4)

Product Code	CSB-BP873654HU
Abbreviation	SERTAD4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NUC0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTLVLSMNRF CEPIVSEGAA EIAGYQTLWE ADSYGGPSPP GPAQAPLQGD RGAGPPLAGS HYRGISNPIT TSKITYFKRK YVEEEDFHPP LSSCSHKTIS IFEERAHILY MSLEKLKFD DPEVYLRRSV LINNLMKRIH GEIMQNNWC FPACSFNGTS AQEWFMAQDC PYRKRPMAK EECEKFHACC FYQECGGHYL NLPLSVNANV GSASTAASSP SASSSSSSSS SSPPLPLPSC SRQVDFDVGS ASIYKSDGQI PANEIFVTNV RSLGVQEKAK LNDEKANDDT NRDGGPLSHE PVGNDLAFEC KGQFYDYFET GYNERNNVNE SWKKSLRKKE ASPPSNKLCC SKGSKI
Source	Baculovirus
Target Names	SERTAD4
Protein Names	Recommended name: SERTA domain-containing protein 4
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.