



Recombinant Human CYFIP-related Rac1 interactor B (CYRIB)

Product Code	CSB-YP873658HU
Abbreviation	Recombinant Human CYRIB protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NUQ9
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGNLLKVLCTDLEQGPNFFLDFENAQPTESEKEIYNQVNVVLKDAEGILEDLQ SYRGAGHEIREAIIQHPADEKLEKAWGAVVPLVGLKFKFYEFSSQRLEAALRGL LGALTSTPYSPTQHLEREQALAKQFAEILHFTLRFDELKMTNPAIQNDFSYYRR TLSRMRINNVPAEGENEVNNELANRMSLFYAEATPMLKTLSDATTKFVSENKN LPIENTTDCLSTMASVCRVMLETPEYRSRFTNEETVSFCLRVMVGVIIYDHSVH PVGAFAKTSKIDMKGCIVLKDQPPNSVEGLLNALRYTTKHLNDETTSKQIKSM LQ
Source	Yeast
Target Names	CYRIB
Protein Names	Recommended name: Protein FAM49B Alternative name(s): L1
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.