



Recombinant Human WD repeat-containing protein 70 (WDR70)

Product Code	CSB-MP873668HU
Abbreviation	WDR70
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NW82
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MERSGPSEVT GSDASGPDPQ LAVTMGFTGF GKKARTFDLE AMFEQTRRTA VERSRKTLEA REKEEEMNRE KELRRQNEDI EPTSSRSNVV RDCSKSSSRD TSSSESEQSS DSSDDELIGP PLPPKMVGKP VNFMEEDILG PLPPPLNEEE EEAEAAAAAAAA EEEENPVHKI PDSHEITLKH GTKTVSALGL DPSGARLV TG GYDYDVKFWD FAGMDASFKA FRSLQPCECH QIKSLQYSNT GDMILVVS GS SQAKVIDRDG FEVMCEIKGD QYIVDMANTK GHTAMLHTGS WHPKIKGEFM TCSNDATVRT WEVENPKKQK SVFKPRTMQG KKVIPPTCTY SRDGNLIAAA CQNGSIQIWD RNLTVHPKFH YKQAHDSGTD TSCVTFSYDG NVLASRGGDD SLKLWDIRQF NKPLFSASGL PTMFPMTDCC FSPDDKLIVT GTSIQRGCGS GKLVFFERRT FQRVYEIDIT DASVVRCLWH PKLNQIMVGT GNGLAKVYYD PNKSQRGAKL CVVKTQRKAK QAETLTQDYI ITPHALPMFR EPRQRSTRKQ LEKDRLDPLK SHKPEPPVAG PGRGGRVGTH GGTLSYIVK NIALDKTDDS NPREAILRHA KAAEDSPYWV SPAYSKTQPK TMFAQVESDD EEAKNEPEWK KRKI
Source	Mammalian cell
Target Names	WDR70
Protein Names	Recommended name: WD repeat-containing protein 70
Expression Region	1-654
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.