



Recombinant Human BTB/POZ domain-containing protein KCTD5 (KCTD5)

Product Code	CSB-BP873685HU
Abbreviation	KCTD5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NXV2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AENHCELLS PARGGIGAGL GGGLCRRCSA GLGALAQRPG SVSKWVRLNV GGTYFLTTRQ TLCRDPKSFL YRLCQADPDL DSDKDETGAY LIDRDPTYFG PVLNLYLRHGK LVINKDLAEE GVLEEAEFYN ITSLIKLVKD KIRERDSKTS QVPVKHVYRV LQCQEEELTQ MVSTMSDGWK FEQLVLSIGSS YNYGNEDQAE FLCVVSKELH NTPYGTASEP SEKAKILQER GSRM
Source	Baculovirus
Target Names	KCTD5
Protein Names	Recommended name: BTB/POZ domain-containing protein KCTD5
Expression Region	2-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.