



Recombinant Human EH domain-containing protein 2 (EHD2)

Product Code	CSB-MP873710HU
Abbreviation	EHD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZN4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFSWLKRGG A RGQQPEAIR T VTSALKELYR TKLLPLEEHY RFGAFHSPAL EDADFDGKPM VLVAGQYSTG KTSFIQYLLE QEVPGSRVGP EPTTDCFVAV MHGDTEGTVP GNALVVDPDK PFRKLNPFNG TFLNRFMCAQ LPNQVLESIS IIDTPGILSG AKQRVSRGYD FPAVLRWFAE RVDLIILLFD AHKLEISDEF SEAIGALRGH EDKIRVVLNK ADMVETQQLM RYVGALMWAL GKVVGTPEVL RVYIGSFWSQ PLLVPDNRRL FELEEQDLFR DIQGLPRHAA LRKLNDLVKR ARLVRVHAYI ISYLKKEMPS VFGKENKKKQ LILKLPVIFA KIQLEHHISP GDFPDCQKMQ ELLMAHDFTK FHSLKPKLLE ALDEMLTHDI AKLMPLLRQE ELESTEVGVQ GGAFEGTHMG PFVERGPDEA MEDGEEGSDD EAEWVVTCKD SKYDEIFYNL APADGKLSGS KAKTWMVGTK LPNSVLGRIW KLSDVDRDGM LDDEEFALAS HLIEAKLEGH GLPANLPRL VPPSKRRHKG SAE
Source	Mammalian cell
Target Names	EHD2
Protein Names	Recommended name: EH domain-containing protein 2 Alternative name(s): PAST homolog 2
Expression Region	1-543
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the EH domain-containing protein family. These proteins are characterized by a C-terminal EF-hand domain, a nucleotide-binding consensus site at the N terminus and a bipartite nuclear localization signal. The encoded protein interacts with the actin cytoskeleton through an N-terminal domain and also binds to an EH domain-binding protein through the C-terminal EH domain. This interaction appears to connect clathrin-



dependent endocytosis to actin, suggesting that this gene product participates in the endocytic pathway.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.