



Recombinant Alligator mississippiensis L-lactate dehydrogenase B chain (LDHB)

Product Code	CSB-YP874049AZO
Abbreviation	LDHB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9PW05
Product Type	Recombinant Protein
Immunogen Species	Alligator mississippiensis (American alligator)
Purity	>85% (SDS-PAGE)
Sequence	STLKEKLIT QIASERTIPN SKVTVVGVGQ VGMACAISIL GKRLGDELAL VDVWEDKLLKG EMMDLQHGS LFLQTHKIVAD KDYAVTANSK IVVVTAGVRQ QEGESRLNLV QRNVNVFKFI IPQIIKYSPD CTILVVSNPV DILTYVTWKL SGLPKHRVIG SGCNLD SARF RYLMSEKLG IHPSSCHGWIL GEHGDSSVAV WSGVNVAGVS LQELNPAMGT DRDSEKWKEV HKQVVESAYE VIKLLKGYTNW AIGLSVADLL ETMMKNLCRV HPVSTMVKGM YGIENEVFLS LPCVLSASGL TSVINQKLKD DEVAQLRSSA DTLWSIQKDL KDL
Source	Yeast
Target Names	LDHB
Protein Names	Recommended name: L-lactate dehydrogenase B chain Short name= LDH-B EC= 1.1.1.27
Expression Region	2-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.