



Recombinant Human immunodeficiency virus type 1 group M subtype K Protein Nef (nef)

Product Code	CSB-BP874074HLQ
Abbreviation	nef
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9QBZ7
Product Type	Recombinant Protein
Immunogen Species	Human immunodeficiency virus type 1 group M subtype K (isolate 97ZR-EQTB11) (HIV-1)
Purity	>85% (SDS-PAGE)
Sequence	GGKWSKSSI VGWSTVRERM RKTTPAADGV GAVSQDLDKH GAVTSSNTAF NNPDCAWLEA QEDEDVGFPV RPQVPLRPMT FKGAFDLGFF LKEKGGLDGL IYSKRRQEIL DLWVYHTQGF FPDWQNYTPG PGIRYPLTFG WCYKLVVPDP REVEEATEGE NNCLLHPVNQ HGMEDEHREV LKWKFDSSLA RKHVAREMHP EYYKDC
Source	Baculovirus
Target Names	nef
Protein Names	Recommended name: Protein Nef Alternative name(s): 3'ORF Negative factor Short name= F-protein Cleaved into the following chain: 1. C-terminal core protein
Expression Region	2-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.