



Recombinant *Zymomonas mobilis* Lipoyl synthase (lipA)

Product Code	CSB-YP874243ZBM
Abbreviation	lipA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9RNY7
Product Type	Recombinant Protein
Immunogen Species	<i>Zymomonas mobilis</i> subsp. <i>mobilis</i> (strain ATCC 31821 / ZM4 / CP4)
Purity	>85% (SDS-PAGE)
Sequence	MSKIPAPDEK TSRSDSGSSI LRKPDWLRVR SPGGAAFNET HGLIRKLGLA TVCEEAACPN IGECWTKKQA TVMILGEVCT RACAFNCVKT GHPDKVNP PGHVADVAE MNLEHIVITS VDRDDLEDGG ASQFVKVIEA VRARTPKTTI EILTPDFRGK PDHALDMIVK ARPDVFNHNL ETVPRLYPTI RPGARYFTSL RLLEKVRERD PSIFTKSGLM LGLGEDRLEV HQVMDDMRQA DVDFLTGLGQY LQPTPRHVRV EEFVTPDSFK AYAATARAKG FSMVASSPLT RSSYYAGADF ARLKEARQKK LEKIAAKNSG KGA
Source	Yeast
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.