



# Recombinant Arabidopsis thaliana Aspartic proteinase-like protein 2 (At1g65240)

<b>Product Code</b>	CSB-EP874357DOA-B
<b>Abbreviation</b>	At1g65240
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9S9K4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	NFVFKAQ HKFAGKKKLN EHFKSHDTRR HSRMLASIDL PLGGDSRVDS VGLYFTKIKL GSPPKEYHVQ VDTGSDILWI NCKPCPKCPT KTNLNFRLSL FDMNASSTSK KVGCDDDFCS FISQSDSCQP ALGCSYHIVY ADESTSDGKF IRDMLTLEQV TGD LKTGPLG QEVVFGCGSD QSGQLGNGDS AVDGVMGFGQ SNTSVLSQLA ATGDAKRVFS HCLDNVKG GGG IFAVGVDSP KVKTTMPVPN QMHYVMLMG MDVDGTSLDL PRSIVRNGGT IVDSGTTLAY FPKVLVDSL I ETILARQPVK LHIVEETFQC FSFSTNVDEA FPPVSFEFED SVKLTVPYPHD YLFTLEEELY CFGWQAGGLT TDERSEVILL GDLVLSNKL V VYDL DNEVIG WADHNCSSSI KIKDGSGGVY SVGADNLSS
<b>Source</b>	E.coli
<b>Target Names</b>	At1g65240
<b>Protein Names</b>	Recommended name: Aspartic proteinase-like protein 2 EC= 3.4.23.-
<b>Expression Region</b>	24-449
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.